

Site: Dallas
LC/MS/MS

VLX2 AND SCIEX LC/MS/MS MAINTENANCE LOG

Form ID: DALTXF110[1.0]

Instrument:

Bumblebee

Month/Year

"*" =NOT IN USE "/" =NOT RECORDED "X" CHECKED/PERFORMED

Daily Maintenance	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Reboot system PC																															
Check Mobile Phase volumes																															
Check Wash Solvent volumes																															
Check Waste Container																															
VLX2 - "System All On"																															
Equilibrate columns w/ method 20 mins																															
Equil MS w/method 10 mins																															
Perform Column Test Run and verify acceptability																															
Column Pressure (bar)	S1																														
	S2																														
Injection Count	S1																														
	S2																														
Check for leaks																															
Tech Initial																															
As Needed	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Change MPA/MPB ("Change All Bottles")																															
Purge pumps																															
Clean curtain plate																															
Change Needle Wash																															
Change Rear Seal Wash																															
Perform CAM track maintenance																															
Change stainless steel kk	S1																														
	S2																														
Change Column (Note batch	S1																														
and serial number for each stream in comments)	S2																														
Clean Q Jet																															
Clean Orifice																															
Replace Electrode																															
Clean Q0																															
Tech Initial																															

Monthly Review:

Other as Needed Actions Performed:

Comments:

Monthly Review: